

EyeSuite Holladay Setup Quick Guide

Click on the Windows Icon located in the bottom Left Corner of the Desktop.

Click on "Computer" then double click "OS (C:)"

Right Click and Select "New ➡ Folder"

Label the Folder "Holladay IOL Consultant Data"

There may be a folder already labeled "HolladayII" or "Holladay II"

If this is the case then do not worry about creating a new folder as EyeSuite will connect to the existing folder/location.

From the Haag-Streit EyeSuite Main Screen Select "Tools ➡ Preferences ➡ +Biometry ➡ Exports"

Make sure the mode is set to "Automatic"

From the Dropdown menu next to "Format" make sure "Holladay IOL Consultant Connector" is selected

The "Template" and "Filename" will autopopulate.

For "Output Folder" browse to "C:\Holladay IOL Consultant Connector" or the existing "*HolladayII* or *Holladay II*" folder

Then Press the "OK" button at the bottom right corner.

